



WHOLE EXOME SEQUENCING

- ▶ Whole Exome Sequencing (WES) is a comprehensive genetic test that identifies DNA changes in the entire protein-coding region of the genome.
- ▶ WES may contribute to diagnosis of rare diseases in less time and at a lower overall cost than traditional DNA testing.
- ▶ WES is recommended
 - when symptoms are complex and unspecific
 - when a diagnosis of significant molecular heterogeneity is indicated
 - in case of neurodevelopmental disorders, including intellectual disability, global developmental delay, and autism spectrum disorder.



RARE DISEASES

- ▶ Rare diseases are defined in Europe as those that have a prevalence of one in 2000 or less.
- ▶ Over 8000 rare diseases are known affecting approximately 6% of the global population.
- ▶ 70% of rare diseases are genetically determined, the majority of which may be recognized with WES.
- ▶ Without timely, comprehensive genetic testing, patients often have to wait years for the correct diagnosis often following a stressful diagnostic odyssey.



METHOD

- ▶ Isolated genomic DNA or peripheral blood sample in EDTA.
- ▶ Sequencing of total coding regions as well as 3'/5' UTRs of the entire exon.
- ▶ Bioinformatic identification of single nucleotide variants (SNV), short insertions and deletions.
- ▶ Variant classification and annotation (e.g. ClinVar, COSMIC, HGMD, dbSNP, VarSome).


COVERAGE: >95% (>500×)

AVERAGE SEQUENCING DEPTH: >1000×

TURNAROUND TIME: 2-3 WEEKS

WORKFLOW

Our mission is advancing scientific research in the fields of **BIOTECHNOLOGY** and **MEDICINE** as well as applying the latest innovative technologies in diagnostics. **IBIOSCIENCE LTD.** in collaboration with **UNIVERSITY OF PÉCS SZENTÁGÓTHAI RESEARCH CENTER** provides state-of-the-art next generation sequencing services and expertise for the Hungarian scientific community.



REFERENCES

Nambot S et al. Genet Med. 2018 Jun;20(6):645-654. PMID: 29095811
 Helbig KL et al. Genet Med. 2016 Sep;18(9):898-905. PMID: 26795593
 Friedman JM et al. JAMA. doi: 10.1001/jama.2020.11126. PMID: 32735311
 Nguengang et al. Eur J Hum Genet. 2020 Feb;28(2):165-173. PMID: 31527858
 Wagner M et al. J Inherit Metab Dis. 2019 Sep;42(5):909-917. PMID: 31059585
 Lalonde E et al. Front Pediatr. 2020 Jul 8;8:373. PMID: 32733828



CONTACT

order@ibioscience.hu
+36 70 674 6611



SAMPLE DISPATCH

DNA or peripheral blood



SAMPLE PROCESSING

paired-end 150 bp sequencing



GENETIC ANALYSIS

bioinformatic analysis, variant identification and annotation



REPORT

categorization of variants based on guidelines (pathogenic, likely pathogenic, VUS etc.)

iBioScience Ltd.

Dr. Majorossy Imre Str 36., Pécs, 7625, Hungary
 Phone: +36 70 674 6611
 E-mail: order@ibioscience.hu

ibioscience.hu